## labmodul



# SYSTEMACITY

Spacious with extra storage capacity

Prepared for trays and baskets

High quality for daily use

Designed for GMO facilities

Pure design

Great carrying capacity for heavy equipment

Easy to clean





#### System Capacity

CAPACITY fully exploits the room beneath the worktop and offers a large storage capacity.

It is stationary, stable and particularly suited for laboratories with heavy equipment.


The cabinets comes in two heights adjusted to standing or sitting work tasks.

The cabinets are prepared for mounting of 60/40 storage system panels and can be delivered with system trays and baskets – either

Dishwashers and other equipment up to a height of 87 cm can be placed under the bench tops.

CAPACITY can be combined with our other furniture systems.

Benches can be fitted with a service bridge containing light and builtin installations for power, gasses andwater etc.





# D-VG-ABB HHSD-VG2-ABS HSD-RLAB

#### A wide programe

Capacity is offered with cabinet box and fronts made of compact laminate for working environments demanding a solution of inorganic materials, resistant to humidity and water.

Base cabinets are mounted on a laminated waterproof plinth. This is an easy-to-clean solution meeting all GMO furnishing demands.

The cabinets are prepared for the extensive ISO tray/basket system 60/40.

The flexibility provides extra storage room and the tray system offers unlimited possibilities for further division and organization of small parts.

HSS-HG2





#### Capacity for GM

Capacity is suitable for laboratories with high cleaning demands.

Undercabinets are mounted on a waterproof plinth clad in laminate.

This solution covers all GMO requirements.

The cabinets are prepared to adapt to the comprehensive ISO tray/basket system 60/40.

#### Solutions for Clean Rooms

Capacity is offered with cabinet body and -fronts made from compact laminate for clean rooms and other work environments that requires unorganic materials that can withstand water and moisture.





#### Service Bridges

The service bridges/panels have built-in working light that illuminates the worktops with a convincing 900 lux. The distribution of the light above the worktops creates a clear and comfortably lit working area.

The lighting fittings have tightfitting transparent lids for easy and efficient cleaning.

The service bridges/panels are offered in two versions- a low model with horizontal plugs and a high model with vertical plugs.

The low model maximises between opposite working areas.

The high model gives easy access to the plugs and is advantageous where equipment is often plugged in or out.



#### Installations


Both models are prepared for additional mounting of extra outlets.

The service bridges/panels are delivered with all technical installations and outlets for power and plumbing ready for use.

The installations end in a chosen pick-up point for quick and easy connection to the building's installations.



By utilizing the combinations units between 1 and 2 module cover plates, in conjunction with the movable cover can technique tower simply and flexibly meet a variety of functional requirements and configurations.



cover sheets removability gives easy access during installation and servicing. Future modifications / changes are made

easy.



#### **Fronts**


- Heighpressure laminate
- Compact laminate
- Linolium
- Corian



#### Glass cabinet fronts

- Natural anodised 20 mm aluminum frame with 110° hinges
- White powder coated 20 mm aluminum frame with 110° hinges
- Laminate 60 mm frame with 270° hinges
- Laminate 60 mm frame concealed 110° hinges





#### Handles

- Standard logistics aluminum handles with inclined surface for for example Dymo strips.
   Width 190 mm
- Handle in brushed stainless steel.
  Width 160 mm
- Handle in matt aluminium.
  Width 160 mm







#### Countertops

#### **Material options**

- Laminate
- Compact laminate
- Stainless steel
- Corian
- Glass

#### Laminate colors





#### **Fronts**

#### Materials

- Laminate
- Compact laminate
- Corian

#### Laminate colors



#### **Plinths**

Two standard colors







#### Tabletop profiles



Rear edge

Postform





Compact laminate with Alu upturned edge





Trespa + laminate with Alu upturned edge





Stainless steel, with bead





Epoxy Resin with bead





Polypropylen with bead









#### Accessories for Service Bridges

The unique bridge is prepared for various mounting solutions that intelligently indemnify the desktop for other tasks. The suspensions can be moved throughout the length of the bridge so that the desired function is at your fingertips



TOWEL HOLDER Aluminium finish, width 570 mm

MOUNTINGS STAND Aluminium finish, width 380 mm



PAPERTOWELHOLDER Aluminium finish, W352 x D150 x H255 mm

PAPERTOWELHOLDER - 3 levels Aluminium finish, W352 x D150 x H305 mm

KNIFE HOLDER Aluminium finish, W380 x D80 x H370 mm



UNIVERSAL SHELF Aluminium finish, W440 x D150 x H275 mm

DOUBLE SHELF Aluminium finish, W380xD80xH275 mm

DOCUMENT HOLDER Aluminium finish, W440 x D74 x H320 mm




DEVIDER HOLDER Aluminium finish, W380 x D190 x H255 mm

WASH RACK Aluminium finish, W380 x D80 x H370 mm

CUP HOLDER Aluminium finish, W147 x D150 x H335 mm, Incl. cup  $\emptyset$ 147 mm

#### Sink worktops



#### Polypropylene sinks

PPDK1	265x117x185 mm, Drip cup
PPDK2	Ø170x115 mm, Drip cup
PP320	320x320x200 mm
PP400	400x400x250 mm
PP500	500x400x250 mm
PP600	600x400x250 mm

#### Worktop materials

- High-pressure laminate
- Compact laminate
- Polypropylene (PP)
- Epoxy Resin
- Acid resistant stainless steel ASI316
- Ceramics




#### Ceramic sinks

KKDK1	Ø105x112 mm, Drip cup
KKDK2	250x95x112 mm, Drip cup
KK400	400x400x250 mm
KK500	500x400x250 mm
KK600	600x400x250 mm






Drain plug 97B0070P



Overflow pipe 158mm 97B0158





**Epoxy Resin sinks** 


\$497 150x70x108 mm, Drip cup \$499 229x76x134 mm, Drip cup 250x125x150 mm, Drip cup \$1000 SDS5-DI 300x200x200 mm 356x254x152 mm SA5-DI SN9-DI 406x305x203 mm SGL 10-DI 400x400x300 mm SM8-DI 450x365x200 mm SEL 16-DI 600x450x300 mm SVL15-DI 635x381x254 mm





TR100	Ø100 mm round Funnel
TR150	Ø150 mm round Funnel
DK130	Ø130mm round Drip cup
DK250	250x150x100 mm, Drip cup
A340	340x300x200 mm
A400	400x340x145 mm
A480	480x340x180 mm

A500 500x400x160 mm A600 600x400x200 mm

#### Production information








#### Cabinets

DEPTH

600


800


Top cabinet - Wall hanging




WIDTH

600






WIDTH



WIDTH

900/1200

Under bench cabinet on plinth - low






8 8




#### Tall storage cabinets

#### Storage cabinets HSS




#### Adjustable shelves







#### Storage cabinets system 60/40 for trays and baskets



Trays and baskets




Storage cabinets KSS for chemicals and poison with moveable shelfs protected by glass / depth 450 mm and 600 mm




Adjustable shelves with protective glass






Storage cabinets KSS for chemicals and poison with pull-out / depth 600 mm




Pull-out door and 5 pc. of moveable wire baskets



Pull-out door and drawer with soft-close




Fire cabinet / depth 617 mm





Fire Resistive acid / chemicals cupboard with self closing door














### labmodul