

SYSTEM STANDARD

Flexible and movable
Future-orientated and suited for all purposes
High quality for daily use
Great load bearing
Smart design
High vertical
Easy to Clean




SYSTEM STANDARD

STANDARD is a movable furniture system with strong constructions designed for easy combination. The system is compatible to our other system solutions which.

All constructions are thoroughly tested and the furniture system is produced of resistant and robust materials to encure a durable solution for many years of daily use. Furthermore, the materials are carefully chosen with focus on recycling and sustainability.

The STANDARD furniture system is exceedingly flexible because of the aluminium framework that makes it possible to shorten the modules to the desired width.



A broad program


Standard is delivered with a broad range of different worktop and front materials to choose from along with different types of hinges and drawer handles to fit the desired cabinet function precisely.

Drawers are delivered with 80 % soft close pull out as standard or with 110 % soft close + pull out.

The base cabinets are mounted hanging in runners and can easily be moved by pulling the cabinet sideways

The flexibility of the system allows seating, workstations and sub worktable equipment to be established where needed.









Good ergonomics

The tables can be delivered with different types of electrical heightadjustability.

We offer 2 different control panels for the height-adjustability systems: up / down buttons

up/down buttons with memory function for 3 personal settings







Activator for normal loads (200 kg.)



Hydraulic for heavy loads (300 kg.)



Mobilekit loads (150 kg.)





Service Bridges


The service bridges/panels have built-in working light that illuminates the worktops with a convincing 900 lux. The distribution of the light above the worktops creates a clear and comfortably lit working area. The lighting fittings have tight-fitting transparent lids for easy and efficient cleaning. The service bridges/panels are offered in two versions- a low model with horizontal plugs and a high model with vertical plugs. The low model maximises between opposite working areas. The high model gives easy access to the plugs and is advantageous where equipment is often plugged in or out.

Installations


Both models are prepared for additional mounting of extra outlets.

The service bridges/panels are delivered with all technical installations and outlets for power and plumbing ready for use.

The installations end in a chosen pick-up point for quick and easy connection to the building's installations.



By utilizing the combinations units between 1 and 2 module cover plates, in conjunction with the movable cover can technique tower simply and flexibly meet a variety of functional requirements and configurations.



Cover sheets removability gives easy access during installation and servicing. Future modifications / changes are made easy.



Fronts

- Heighpressure laminate
- Compact laminate
- Linolium
- Corian



Glass cabinet fronts

- Natural anodised 20 mm aluminum frame with 110° hinges
- White powder coated 20 mm aluminum frame with 110° hinges
- Laminate 60 mm frame with 270° hinges
- Laminate 60 mm frame concealed 110° hinges





Handles

- Standard logistics aluminum handles with inclined surface for for example Dymo strips. Width 190 mm
- Handle in brushed stainless steel.
 Width 160 mm
- Handle in matt aluminium.
 Width 160 mm







Countertops

Material options

- Laminate
- Compact laminate
- Stainless steel
- Corian
- Glass

Fronts



Corian

Laminate colors



Laminate colors

Laminate colors


White



Plinths

Two standard colors

Tabletop profiles



Terril

Polypropylen with bead

Ceramic with bead



Accessories for Service Bridges

The unique bridge is prepared for various mounting solutions that intelligently indemnify the desktop for other tasks. The suspensions can be moved throughout the length of the bridge so that the desired function is at your fingertips



TOWEL HOLDER Aluminium finish, width 570 mm

MOUNTINGS STAND Aluminium finish, width 380 mm



PAPERTOWELHOLDER Aluminium finish, W352 x D150 x H255 mm

- PAPERTOWELHOLDER 3 levels Aluminium finish, W352 x D150 x H305 mm
- KNIFE HOLDER Aluminium finish, W380 x D80 x H370 mm



UNIVERSAL SHELF Aluminium finish, W440 x D150 x H275 mm



DEVIDER HOLDER Aluminium finish, W380 x D190 x H255 mm



DOUBLE SHELF Aluminium finish, W380xD80xH275 mm

1

WASH RACK Aluminium finish, W380 x D80 x H370 mm



DOCUMENT HOLDER Aluminium finish, W440 x D74 x H320 mm



CUP HOLDER Aluminium finish, W147 x D150 x H335 mm, Incl. cup Ø147 mm

Sink worktops



Worktop materials

- High-pressure laminateCompact laminate
- Polypropylene (PP)
- Epoxy Resin
- Acid resistant stainless steel ASI316
- Ceramics



156



Overflow pipe 158mm 97B0158



Water trap 97B0080




Polypropylene sinks

PPDK1 265x117x185 mm, Drip cup PPDK2 Ø170x115 mm, Drip cup PP320 320x320x200 mm 400x400x250 mm PP400 PP500 500x400x250 mm PP600 600x400x250 mm

Ceramic sinks

KKDK1	Ø105x112 mm, Drip cup
KKDK2	250x95x112 mm, Drip cup
KK400	400x400x250 mm
KK500	500x400x250 mm
KK600	600x400x250 mm



Epoxy Resin sinks

S497	150x70x108 mm, Drip cup
S499	229x76x134 mm, Drip cup
S1000	250x125x150 mm, Drip cup
SDS5-DI	300x200x200 mm
SA5-DI	356x254x152 mm
SN9-DI	406x305x203 mm
SGL 10-D	01 400x400x300 mm
SM8-DI	450x365x200 mm
SEL 16-D	600x450x300 mm
SVL15-DI	635x381x254 mm




Stainless steel sinks

A340 340x300x200 mm A400 400x340x145 mm A480 480x340x180 mm A500 500x400x160 mm A600 600x400x200 mm	A400 A480 A500	400x340x145 mm 480x340x180 mm 500x400x160 mm
---	----------------------	--

labmodul

Product Information



labmodul

USO-D-001

USO-V-001

USO-H-001

USO-SA-031

USO-SB-033


USO-SC-035

USO-SD-037




Tall storage cabinets


Storage cabinets HSS



Storage cabinets system 60/40 for trays and baskets









2000

100







WIDTH

Trays and baskets








See the selection of trays and baskets in a separate brochure


labmodul




Storage cabinets KSS for chemicals and poison with moveable shelfs protected by glass / depth 450 mm and 600 mm

Adjustable shelves with protective glass






Storage cabinets KSS for chemicals and poison with pull-out / depth 600 mm






Pull-out door and 5 pc. of moveable wire baskets






Pull-out door and drawers with soft-close








Fire Resistive acid / chemicals cupboard with self closing door



labmodul


Product Information

Wall mounted





Service Bridges


Freestanding sections



Installations









Valhøjs Alle 190, DK-2610 Rødovre, Tlf. +45 70 21 21 89